



Clinivision 49GT 4K

Article	9200810
Determination and possibility of application	Display unit with optional computer module for use in hygienically sensitive areas, viewing of medical images/video recordings, digital radiological examinations, documentation in HIS/RIS/PACS or PDMS
Mounting Options	Inwall mounted and surface-mounted mounting, mobile solution
Testing / Conformity	EU Regulation 2017/745, MPG Class I
Warranty	5 years according to our warranty policy

Display

Type Display	IPS TFT LED
Display Size	49" (1232 mm diagonal)
Resolution	3840 x 2160 pxl (4K)
Color depth	10bit/1.07 billion Colours
Brightness / Contrast	500 cd/m ² / 1100:1
Perspective	R/L 178, O/U 178
Reaction time	8 ms Grey to Grey
Signaleingänge	1 x DVI-D, 1 x Display Port, Optionen: PiP, loop through
Brightness stabilisation	ALS (<i>ambient light sensor</i>) brightness adjustment optional
OSD	OSD intern, IR-Control optional
LUTE	DICOM, GAMMA 1.8/2.2, user specified LUT

Case

Protection	IP65 Front
Basic housing	Anodized aluminum housing
Frontglas	5 mm ESG anti reflective, Floatglas
Touch	PCAP optional
Keyboard Tray	infinitely adjustable keyboard/mouse tray with integrated USB, optional clinitouch keyboard and mouse
Device dimensions	B 1299 x H 979 x T 153,40 mm
Weight	approx. 59 kg

General

Power supply	internal power supply, AC 90-260V (47-63Hz)
Electricity consumption	98 VA (without computer module)



Clinivision 49GT 4K

