



clinvision 24GT

Article	9200982
Determination and possibility of application	Display unit with optional computer module for use in hygienically sensitive areas, viewing of medical images/video recordings, digital radiological examinations, documentation in HIS/RIS/PACS or PDMS
Montagemöglichkeit	Surface-mounted and inwall mounting
Testing/Conformity	EU Regulation 2017/745, MPG Class I
Warranty	5 years according to our warranty policy

Display

Type Display	TFT LED
Display Size	60,5 cm (24")
Resolution	1920 X 1080/16:09
Pixel Pitch	0.274 mm x 0.274 mm
Color depth	10bit / 1.07 billion Colours
Brightness / Contrast	500 cd/m ²
Perspective	R/L 178, O/U 178
Reaction time	14 ms Grey to Grey
Signal inputs	1 x DVI-D, 1 x Display Port (for PC module)
Brightness Stabilisation	ALS (<i>ambient light sensor</i>) brightness adjustment optional
OSD	OSD intern, IR-Control optional
LUTE	DICOM, GAMMA 1.8/2.2, user specific LUT

Case

Protection	IP65 Front
Basic housing	Anodized aluminium housing
Frontglas	4,0 mm ESG antireflective Floatglas
Touch	PCAP optional
Input System	Infinitely adjustable keyboard/mouse tray with integrated USB, optional clintouch keyboard and mouse
Device dimensions	min. B700xH730xT100cm
Weight	approx. 31 kg

General

Power supply	Internal power supply, AC 90-260V (47-63Hz)
Electricity consumption	max 28 VA