



clinvision 55GT 4K

Article	9200982
Determination and possibility of application	Display unit with optional computer module for use in hygienically sensitive areas, viewing of medical images/video recordings, digital radiological examinations, documentation in HIS/RIS/PACS or PDMS
Montagemöglichkeit	Surface-mounted and inwall mounting
Testing/Conformity	EU Regulation 2017/745, MPG Class I
Guarantee	5 years according to our warranty policy

Display

Type Display	IPS TFT LED
Display Size	55" (1397 mm diagonal)
Resolution	3840 x 2160 pxl (4K)
Color depth	10bit / 1.07 billion Colours
Brightness / Contrast	500 cd/m ²
Perspective	R/L 178, O/U 178
Reaction time	8 ms Grey to Grey
Signaleingänge	1 x DVI-D, 1 x Display Port
Helligkeitsstabilisation	ALS (<i>ambient light sensor</i>) brightness adjustment optional
OSD	OSD intern, IR-Control optional
LUTE	DICOM, GAMMA 1.8/2.2, user specific LUT

Case

Protection	IP65 Front
Basic housing	Anodized aluminum housing
Frontglas	6,0 mm ESG antireflective Floatglas
Touch	PCAP optional
Eingabesystem	Infinitely adjustable keyboard/mouse tray with integrated USB, optional clintouch keyboard and mouse
Device dimensions	B 1390 x H 1030 x 155 mm
Weight	approx. 89 kg

General

Power supply	Internal power supply, AC 90-260V (47-63Hz)
Electricity consumption	max 95 VA